

Need help getting your patients screened for colorectal cancer?

Quest Diagnostics® can ship InSure® ONE™ fecal immunochemical test (FIT) kits directly to your practice. These tests qualitatively detect human hemoglobin from blood in fecal samples, which may be indication of lower gastrointestinal bleeding associated with disorders such as diverticulitis, ulcerative colitis, polyps, colorectal cancers, or large adenomas that bleed.

Is FIT testing accurate?

FIT testing is supported by medical guidelines—including those of the American Cancer Society, the American College of Gastroenterology, and the US Preventive Services Task Force—for early detection of colorectal cancer.

Is the test easy for my patients to use?

InSure ONE is the only FIT test that is performed using toilet water collected from a single bowel movement. There is no stool collection or handling required, and there are no dietary or medicinal restrictions. This helps increase the likelihood your patient will actually complete the test.

What do I need to do?

- Start ordering InSure ONE kits using **Quantum™ Lab Services Manager**. If you don't have an account, create one in less than 5 minutes by clicking [here](#). If you would prefer to order kits over the phone, please call 1.866.MY.QUEST and state “supplies” upon connection for assistance.
- Search for test codes 11290 or 11293 (Medicare) for InSure ONE kits, where you will find all patient instructions for specimen collection and for returning completed test kits.
- Quest will ship the kits to your practice at no cost, based on usage and test completion.
- Distribute kits to average-risk patients aged 50+ during annual wellness visits. A completed requisition form from Quantum Lab Services Manager or your EMR should be printed and stapled to the Reply Form inside the kit before providing to patients.
- You can also view patient-specific test history from the Quantum Lab Services Manager platform.

What does the patient need to do?

- The patient should fill in all items in the patient section on the Reply Form, including patient information, sample information, and insurance billing information. They must also fill out the removable label.
- After reading the collection instructions inside the kit and collecting and placing a water sample on the test card, the completed label should be placed on the card to seal the flap. The test card should then be placed in the mailing envelope provided along with the completed Reply Form.
- Drop into US Mail within 48 hours of collection.

What if I have questions?

If you have any questions or if you need additional assistance, please call **1.866.MY.QUEST** or go to **QuestInSureONE.com** to watch our InSure ONE collection video.